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1 General

1.1 Important Information

Please check immediately after *unpacking* that the contents of the consignment conforms to the information given on the delivery note.

Please read the *operating instructions* before you operate the unit and follow them in all respects.

No liability will be accepted for personal injury nor material damages in the event that damage or breakdowns occur as a result of failure to comply with these operating instructions; neither will any guarantees relating to repairs to or replacements of our products apply.

All units comply with the Federal German Law concerning Technical Implementations. The operating instructions comply with the German Industrial Standard DIN 8418.

If you perform *repair or maintenance work* on units which have come into contact with substances which are detrimental to health, please observe the relevant regulations.

If you return units to us for repair or maintenance work, please follow the instructions below:

- *contaminated units* (radioactively or chemically etc.) must be decontaminated in accordance with the regulations before they are returned.
- Units returned for repair or maintenance must bear a clearly visible note "*Free from harmful substances*". This note must also be provided on the delivery note and accompanying letter.
- Please use the attached attestation declaration.

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"materials and preparations in accordance with the EEC Specification dated 18 September 1979, Article 2"

and in the U.S.A. as

"materials in accordance with the Code of Federal Regulations (CFR) 49 Part 173.240 Definition and Preparation".

Technical modifications reserved.

1.2 Safety Instructions
Diaphragm vacuum pumps represent state-of-the-art technology and are designed to provide optimum operational reliability. The user may, however, be exposed to hazards if pumps are used improperly or for other than their intended purpose.

If pumps are used for any other than their intended purpose all liability and warranty claims will lapse.

Use for the intended purpose also means that the installation, commissioning, operating and maintenance instructions of the manufacturer must be complied with.

- Operating and service personnel must have read and understood the operating instructions for the product.
- All warnings must be observed and all precautions taken.
- The operating and service personnel must be informed of all hazards which might occur in connection with the rotary vane vacuum pump. Personnel must be able to recognize dangers and take preventive measures.
- Any use by unauthorized personnel

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- or careless handling may increase the potential danger.
- The procedures described in the attached operating instructions must be observed in all installation, commissioning, operating and maintenance work.
- No operating modes must be used which may affect the safety of the operator and the pump.
- All unauthorized modifications and alterations affecting safety are prohibited and would cause any liability or guarantee claims to be repudiated.
- Only genuine parts and accessories may be used in service and maintenance work.
- All safety protection measures must be tested on completion of electrical installation and electrical maintenance work (e.g. earthing resistance test).
- All relevant local safety and accident prevention regulations apply for operation of the system.
- Any unclear points with regard to safety, operation and maintenance should be clarified with your nearest Balzers agency or subsidiary company.

1.2.1 Working with Diaphragm Vacuum Pumps; Safety Information

- Pumps must not be used for pumping gases which can cause explosive or inflammable mixtures, nor should they be used where such gaseous mixtures are present.
- Where hot surfaces are involved, allow sufficient space between the position of the pump and such surfaces and protect against the danger of electric sparks.
- The area in which the pump is installed must conform to the official regulations. The pump may only be operated provided protective earthing has been connected.
- All further relevant safety requirements (e.g. accident prevention regulations for laboratories) must be complied with.

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1.3 Application

Diaphragm pumps are noted for their easy servicing, completely oil-free pumping of gases and quiet running.

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Diaphragm pumps are designed for installation in pumping stations. The vacuum connection is via a standard pneumatic hollow bolt onto which can be directly screwed e.g. the fore-vacuum connection of a turbomolecular pump. The exhaust of the standard version is provided with a silencer. Mains connection and thermal motor protection are by means of a one meter long five-core cable without switch and plug. On the 100 V and 110/115 V versions the capacitor is freely suspended on a cable of approx. length 30 cm and must be secured and protected against damage etc. when fitting into the rack.

The "C" version is designed for operations involving corrosive gases.

1.4 Design

The MD 4T is a three stage and the MZ 2T a two stage dry compressing diaphragm vacuum pump which means that the pumped medium does not come into contact with any lubricants.

The pumps are of the displacement type where the movement of the diaphragm causes a periodic change in the size of the working chamber. The flow of gas causes the self-acting valves to open and close.

The pumping units are directly connected to the drive motor.

All pump bearings are arranged outside the working chamber, are enclosed and have lifetime lubrication. Under normal operating conditions the pumps need no servicing.

The MD 4TC Diaphragm Pump is distinguished from the aforementioned standard version by the use of corrosive gas resistant materials in the working chamber which comes into contact with the media (see Section 5.1 for such materials).

To avoid condensation in the pump, the MD 4TC is fitted with a gas ballast valve as standard.

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
Fig. 3.2
MD 4T

Diaphragm pumps A, B, C, D
IN Intake connection
EX Pressure connection

- 1 Intake connection unit
- 2 Silencer
- 3 Exhaust connection
- 28 Hose connection

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2 Installation

2.1 Base and Location

- Position the pump on an even, horizontal surface.
- Fixing to the surface is unnecessary.
- Anchor when fitting in a rack or frame.
- Sight glass 40 must be visible.
- When installing the pump in a closed housing, ensure there is sufficient air circulation.

Caution!

Observe the safety instructions in Section 1.2.1.

2.2 Connection of Pumping Stages

Check the hose connections for the correct type of connecting pieces before connecting the pump to the vacuum chamber.

The intake connections of the pump are marked with "IN" and the outlets with "EX".

Caution !

Interchanging the connections will result in impermissibly high overpressure levels.

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MD 4T, MD 4TC

Fig. 3.1 and 3.2

The pump is equipped with four separately operating diaphragm pumps A, B, C, D which are connected via hoses to form a three-stage vacuum pump.

- The intake connection unit 1 is at A-IN and B-IN.
- The pressure connection from A-EX, B-EX is at intake connection C-IN. This means that A and B operate synchronously. Non-return valves ensure that the two intake and pressure ports are closed to one another so there is no risk of pressure equalisation.
- C-EX is connected to D-IN. Outlet sides 3 of the MD 4TC and the MD 4T are fitted with a DN 10 hose nipple and a silencer respectively.

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MZ 2T

The MZ 2T is a two-stage pump. The intake connection is marked with "IN" and the exhaust connection with "EX" so as to avoid the connections getting mixed up. The outlet side of the MZ 2T is provided with a silencer.

2.3 Intake Connection

- Remove the locking cap from the intake connection.
- Make the connection between pump and vacuum system as short as possible.
- Connect pump to the vacuum system via the intake connection.
- If a vapour releasing liquid is present in the system to be evacuated, a condensate separator must be fitted between pump and pumping system.

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2.4 Exhaust Connection

Important:

Overpressure on the exhaust side must be avoided at all costs (see Section 5, Technical Data for maximum permissible counter pressures at the outlet).

- Lay the exhaust line without shut-off devices.
- If a shut-off device has to be fitted for operating reasons, or if there is a danger that overpressure could build up in the line, official accident prevention regulations must be observed.
- The line from the pump should be laid sloping downwards so that condensate cannot flow back into the pump, otherwise a trap must be fitted.
- Connecting the exhaust side of the pump with a closed system (e.g. the vacuum chamber) must be avoided at all costs otherwise there is a danger the pump could explode.

Caution !

Expelled gases and vapours can represent a health hazard and/or be damaging to the environment.

2.5 Electrical Connections

Diaphragm pumps are driven by alternating current motors.

Important:

Electrical connections must comply with local regulations.

- Ensure that the mains voltage level is compatible with voltage and frequency values information stated on the rating plate.
- The pump may only be connected to mains supply systems which incorporate earthing protection.
- Pumps are equipped with a free embedded thermostatic cut-out in the winding.

Caution !

If the thermostatic cut-out has not been properly wired, the pump has no motor protection.

Diaphragm pumps should be connected in the pumping station as per Wiring Diagram PM 021 557 -S or PM 031 475 -S.

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Legend for Wiring Diagram PM 021 557 -S

- 1 TCW Cooling Water Monitor
- 2 Diaphragm vacuum pump
- 3 Air cooling or TZK Cooling Unit
- 4 Heating
- 8 Cable, TCP turbomolecular pump
- 17 Turbomolecular pump
- 18 Electronic drive unit
- 19 Venting Valve TSF

- A Mains connection
- B Mains connection for remote control

- S1 Pumping station
- S2 Heating
- S4 Pumping station (in remote control)
- S5 Heating (in remote control)
- S6 Reset (in remote control)
- S9 Stand-by
- S10 Stand-by, external

¹⁾ Remove bridge when a TCW is connected

²⁾ 1 - 3 closed when the ambient temperature $\geq 24^{\circ}\text{C}$

³⁾ Remove bridge for remote control operations

⁸⁾ Bridged if a diaphragm vacuum pump is connected without temperature switch

Legend for Wiring Diagram PM 031 478 -S

- 2 Backing pump
- 3 Air cooling
- 18 Electronic drive unit
- 19 Venting valve
- X1 Connection, venting valve
- X2 Connection, air cooling
- X5 Output socket strip TCP
- X8 Mains connection for Electronic Drive Unit TCP 015
- X12 Connection, backing pump
- X25 Connection, electronic drive unit (control signals)
- X26 Connection, backing pump safety valve

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3 Operations

Important:

Before putting the pump into operation it is essential to ensure that impermissibly high pressures cannot arise on the pressure side.

- The pump can be switched on and off at any time.
- Intake and outlet connections are marked with "IN" and "EX" respectively.

Caution !

Interchanging the connections will result in impermissibly high over-pressure levels.

- The pump should be allowed to continue running under atmospheric pressure for several minutes before switching off if it has been subjected to the action of condensates.
- The safety instructions regarding working with diaphragm vacuum pumps contained in Section 1.2.1 continue to apply.

4 Maintenance

Under normal operating conditions the pump requires no maintenance. Pump fluid checks are unnecessary. Valves and diaphragms suffer from the effects of wear and tear. At the latest when the attained pressure values begin to decline, the working chamber, the diaphragms and valves should be cleaned and checked for cracks (see Section 8, Inspection).

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5.1 Medienberührende Werkstoffe

5.1 Substances in Contact with the Media

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6 Servicing Information

6.1 Returning to Service Center for Repair

We request that you contact your local PFEIFFER Service Center should problems arise.

All repair orders are carried out on the basis of our General Terms of Delivery only. For repairs, the warranty period is 6 months. Repair work is performed without checkback if the repair costs do not exceed 50 % of the nominal value of the article.

Please note:

The operating medium in the pump must be drained before you return the pump for repair. When returning for repairs, the important introductory notes contained in these instructions must be observed.

Recommendation:

We recommend that you take out a service contract with your local PFEIFFER Service Center, thereby avoiding possible extended periods of inactivity which could be caused by major damage to the pump.

Should you have any questions, please contact our sales companies or agencies which will be pleased to assist you.

Please see back cover for addresses.

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

7 Troubleshooting


Troubleshooting can be divided into four main categories. The sources of malfunction can be ascribed to a number of possible causes. Malfunction corrective measures depend on the cause.

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

8 Inspection

Fig. 6 to Fig. 8

- Before dismantling the pump, check that the malfunction has not been caused by external factors.
- Whenever performing inspection work, make sure the motor cannot get switched on.
- Where necessary, detach pumps from the system before performing inspection work.
- Only dismantle the pump as far as is necessary in order to correct faults.
- Only use benzine or similar agents for cleaning pump parts. Do not use soluble washing agents.

Important:

Ensure compliance with local regulations when working with solvents.

- Only use genuine parts as per the spare parts list when carrying out inspection work.
- Toxic gases and vapours can arise when carrying out inspection work on corrosive gas version pumps ("C"). Pumps must be decontaminated before inspection. Only qualified personnel are permitted to carry out maintenance and repair work and the regulations must be observed.

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Fig. 6

- 1 Motor
- 2 Silencer
- 3 Cheesehead screw
- 4 Sealing ring
- 5 Bow-shaped handle
- 6 Cheesehead screw
- 7 Housing cover
- 8 Head cover
- 9 Diaphragm spring washer
- 10 Diaphragm support washer
- 11 Housing
- 12 Housing panel
- 13 Excenter bushing with counterweight
- 14 Spacer
- 15 Con-rod
- 16 Double diaphragm
- 17 Inlet and outlet valve seal
- 18 Grooved ball bearing
- 19 Countersunk screw
- 20 Cheesehead screw
- 21 Countersunk screw
- 22 Tolerance ring
- 23 Hose connection (complete)
- 25 Stud with hose nipple

Fig. 6

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Fig. 7

| | |
|-----|-------------------------------------|
| 1 | Housing cover |
| 2 | Head cover |
| 3 | Diaphragm spring washer |
| 4 | Con-rod |
| 5 | Housing |
| 6 | Housing panel |
| 7 | Excenter bushing with counterweight |
| 8 | Cheesehead screw |
| 9 | Diaphragm support washer |
| 10 | Double diaphragm |
| 11 | Inlet and outlet valve seal |
| 12 | Bow-shaped handle |
| 13 | Grooved ball bearing |
| 15 | Silencer |
| 16 | Countersunk screw |
| 17 | Countersunk screw |
| 18 | Cheesehead screw |
| 19 | Tolerance ring |
| 22 | Seal |
| 23 | Motor |
| 24 | Hexagonal screw |
| 25 | Spacer |
| 26 | Cheesehead screw |
| 27 | Connection, intake flange |
| 32 | Housing panel |
| 33a | Cover panel, lower |
| 33b | Cover panel, upper |
| 34 | Raised sheet metal screw |
| 35 | Serrated lock washer |
| 36 | Raised cheesehead screw |
| 37 | Hexagonal nut |
| 38 | Fan |
| 39 | Fan protection |
| 40 | Cover panel |

Fig. 7

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Fig. 8

| | |
|-----|-------------------------------------|
| 1 | Housing cover |
| 2 | Head cover |
| 3 | Diaphragm spring washer |
| 4 | Con-rod |
| 5 | Housing |
| 6 | Housing panel |
| 7 | Excenter bushing with counterweight |
| 8 | Cheesehead screw |
| 9 | Diaphragm support washer |
| 10 | Diaphragm |
| 11 | Inlet and outlet valve seal |
| 12 | Bow-shaped handle |
| 13 | Grooved ball bearing |
| 14 | Connecting screw |
| 17 | Countersunk screw |
| 18 | Cheesehead screw |
| 19 | Tolerance ring |
| 20 | Cylinder bolt |
| 21a | Housing cover (inner part) left |
| 21b | Housing cover (inner part) right |
| 23 | Motor |
| 24 | Hexagonal screw |
| 25 | Spacer |
| 26 | Cheesehead screw |
| 32 | Housing panel |
| 33a | Cover panel, lower |
| 33b | Cover panel, upper |
| 34 | Raised sheet metal screw |
| 35 | Serrated lock washer |
| 36 | Raised cheesehead screw |
| 37 | Hexagonal nut |
| 38 | Fan |
| 39 | Fan protection |
| 40 | Cover panel |

Fig. 8

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8.1 Cleaning and Checking the Working Chamber

8.1.1 MZ 2T

Fig. 6

- Detach the pump from the pumping system and ensure the motor cannot get switched on.
- Screw off hose connection between the pump stages.
- Remove the cheesehead screws 20 (four per pump stage).
- Remove housing cover 7 with head cover 8 taking care to prevent head cover 8 falling out of housing cover 7.
- Check diaphragm pair 16 for damage and replace as necessary (see Section 8.2).
- Remove head cover 8 from housing cover 7 to check valve seals 17.
- Taking care with the position of valve seals 17 (Fig. 9), remove and check for damage, and replacing if necessary.
- Clean all parts.
- Re-assemble pump in reverse order.

8.1.2 MD 4T, MD 4TC

Fig. 7, Fig. 8

- Detach the pump from the vacuum system and ensure the motor cannot get switched on.
- Unscrew connecting hose between the respective pump stages.
- Remove screws 26 and take off housing cover 1.
- Taking care with the position of valve seals 11 (Fig. 9), remove and check for damage, and replacing if necessary.
- Lift up head cover 2 and check diaphragms 10 for wear and tear and replace as necessary (see Section 8.2).
- Clean all parts.
- Re-assemble pump in reverse order. When assembling, the pump should be positioned so that the pump stage to be fitted is on top. The parts are then placed on top into position.

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8.2.1 MZ 2T

Fig. 6

- Dismantle pump stage as per Section 8.1.1.
- Remove screw 19. Caution, this screw has been stuck to make secure against loosening.
- The underlying diaphragm support washer 10 must be removed with the utmost care.

Caution !

It is possible that spacers have been fitted between the drive con-rod and the diaphragm support washer for the purpose of compensation. These spacers must be re-fitted at the same height.

- If it is difficult to part the old diaphragm from support washer 10, apply e.g. naphtha or paraffin.
- If necessary, re-fit spacers under the lower diaphragm support washer 10. Fix in place with a little grease to prevent them falling down.
- Fit new diaphragm 16 and diaphragm support washer 10.

Please note:

The diaphragm is a double diaphragm of total thickness 1.6 mm. The diaphragms are specially selected and are supplied in pairs with a maximum tolerance of 1.6 ± 0.05 mm. Replacement diaphragms must be selected on the same basis.

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- Fit diaphragm support washer 9.
- Place some adhesive (e.g. Loctite 601) under the head of screw 19 and screw in diaphragm and spring washer.
- Further assembly in reverse order.

8.2.2 MD 4T

Fig. 7

- Dismantle pump stages as per Section 8.1.2.
- Remove screw 16. Caution, this screw has been stuck into position and sealed.
- Screw off screws 17 and remove housing panel 6.
- Remove diaphragm 10. The diaphragms are usually stuck to the housing and should be carefully pushed out from the con-rod space by hand.

Caution !

It is possible that spacers have been fitted between the drive con-rod and the diaphragm support washer for the purpose of compensation. These spacers must be re-fitted at the same height. They must not be interchanged with another stage.

- If the old pair of diaphragms 10 cannot be parted from diaphragm support washer 9, apply naphtha or paraffin.
- Clean the parts in naphtha and check for wear and tear, replacing defective parts as necessary.
- If necessary, re-fit spacers under the lower diaphragm support washer 9. Fix in place with a little grease to prevent them falling down.
- Fit diaphragm support washer 9 onto con-rod 4.
- Fit new pair of diaphragms 10.

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Please note:

The diaphragm is a double diaphragm of total thickness 1.6 mm. The diaphragms are specially selected and are supplied in pairs with a maximum tolerance of 1.6 ± 0.05 mm. Replacement diaphragms must be selected on the same basis.

- Fit diaphragm support washer 9.
- Place some adhesive (e.g. Loctite 601) under the head of screw 16 and screw in diaphragm and spring washer.
- Further assembly as per Section 8.1.2.

8.2.2.1 MD 4TC

Fig. 8

- Dismantle pump stages as per Section 8.1.2.
- In addition, the front housing cover panel must be dismantled.
- Using a standard trade face spanner (size 40/4) turn diaphragm support washer 9 thereby loosening diaphragm 10.
- Remove diaphragm spring washer 3 with square connecting screw 14, diaphragm 10 and diaphragm support washer 9.
- Check that square connecting screw 14 is still firmly positioned in diaphragm spring washer 3 and tighten if necessary (left-hand thread).
- Fit new diaphragm 10 between diaphragm spring washer 3 and diaphragm support washer 9.
- Refit diaphragm 10 with spring washer and support washer in reverse order.
- Fit housing cover.
- To fit the PTFE connecting hose, tighten the metal lock nut by hand and then by a half a turn with an Allen key.

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Please note:

When inquiring about or ordering spare parts, please indicate product type and serial number.

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Declaration of Contamination of Vacuum Equipment and Components

The repair and/or service of vacuum components will only be carried out if a correctly completed declaration has been submitted. Non-completion will result in delay.

The manufacturer could refuse to accept any equipment without a declaration.

This declaration can only be completed and signed by authorised and qualified staff:

1. Description of component:

- Equipment type/model: _____
- Code No.: _____
- Serial No.: _____
- Invoice No.: _____
- Delivery Date: _____

2. Reason for return:

3. Equipment condition

- Has the equipment been used?
yes no
- What type of pump oil was used?

- Is the equipment free from potentially harmful substances?
yes (go to section 5)
no (go to section 4)

4. Process related contamination of equipment

- toxic yes no
- corrosive yes no
- microbiological hazard*) yes no
- explosive*) yes no
- radioactive*) yes no
- other harmful substances yes no

*) We will not accept delivery of any equipment that has been radioactively or microbiologically contaminated without written evidence of decontamination!

Please list all substances, gases and by-products which may have come into contact with the equipment:

| Tradename Product name Manufacturer | Chemical name (or Symbol) | Danger class | Precautions associated with substance | Action if spillage or human contact |
|---|------------------------------|--------------|--|--|
| 1. | | | | |
| 2. | | | | |
| 3. | | | | |
| 4. | | | | |
| 5. | | | | |

5. Legally Binding Declaration

I hereby declare that the information supplied on this form is complete and accurate. The despatch of equipment will be in accordance with the appropriate regulations covering Packaging, Transportation and Labelling of Dangerous Substances.

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⇒ **GB, IE**

Declaration of conformity pursuant to the following EU directives:

- Machinery Directive 89/392/EEC
- Electromagnetic Compatibility Directive 89/336/EEC
- Low Voltage Directive 73/23/EEC

We hereby certify that the product specified below is in accordance with the provision of EU Machinery Directive 89/392/EEC - Annex II A, EU Electromagnetic Compatibility Directive 89/336/EEC and EU Low Voltage Directive 73/23/EEC.

The guidelines, harmonised standards, national standards and specifications which have been applied are listed below.

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Guidelines, harmonised standards, national standards in languages and specifications which have been applied:

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Abbreviated Instructions
for
Diaphragm Vacuum Pumps
MZ 2T, MD 4T, MD 4TC

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These abbreviated instructions are only valid together with the detailed operating instructions.

- Detach this page and place near the pump.

Initial Operation

Important:

Only "C" version pumps may be used where operations involve corrosive gases. Comply with the safety instructions.

Installation

- For location, surface requirements and connections please refer to Section 2 of the Operating Instructions.

Caution:

Before operating, remove locking cap and ensure that impermissibly high pressure cannot arise on the pressure side.

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Subsequent Operations

Switching On

- The pump can be switched on at any time. Observe maximum counter-pressure values at the outlet (see Section 5, Technical Data in the operating instructions).
- Under normal operating conditions the pump requires no maintenance.
- If the pump has been subjected to the action of condensates it should be allowed to continue to run at atmospheric pressure for several minutes before switching off.

Connection of the Pump Stages

Check hoses for the correct connections before connecting pumps to the installation.

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Connection must be carried out as per Section 2.2 of the operating instructions.

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Electrical Connections

- Ensure compliance with voltage and frequency information on the rating plate.
- The motor must be connected in compliance with local electrical regulations.

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